

SchureMed Saphenous Vein Harvesting Pad | 508-0385



Stock Photo - Options May Vary

FACTORY NEW

The Saphenous Vein Harvesting Pad (508-0385) positions the patient's legs with external rotation and slight flexion to optimize access for saphenous vein assessment and harvesting during coronary bypass surgery. The pad accommodates minimally invasive, bridging, and open-harvesting procedures.

The pad features a multi-directional stretch cover with internal construction consisting of a conforming foam top layer and firm supportive base. This design provides both comfort and stability during extended procedures.

The water-resistant cover is latex-free with sealed, watertight seams for straightforward cleaning and maintenance between uses.

FEATURES

- Holds legs externally rotated and slightly flexed for minimally invasive, bridging, or open-harvesting procedures
- The cover is a multi-directional stretch cover
- Insides consist of a top layer of conforming foam with a firm, supportive base
- Water-resistant, not made with a natural rubber latex cover
- Sealed, watertight seams for easy clean-up

SPECIFICATIONS

- Pad Dimensions: 33"L x 21"W x 7"H
- Construction: Multi-directional stretch cover over conforming foam and supportive base

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)