

SchureMed Scatter Reducing Surgical Gloves

FACTORY NEW

SchureMed Scatter Reducing Surgical Gloves. These light-touch leaded surgical gloves provide robust protection while allowing maximum flexibility, dexterity, and touch sensitivity. Designed to attenuate both direct and scatter radiation, they offer essential protection during radiographic procedures.

The gloves feature a thin, ergonomic design that enhances freedom of movement, ensuring comfort during extended use. Their superior radiation-attenuating properties make them ideal for professionals requiring both safety and precision. Perfect for surgeons and medical staff working in high-radiation environments, they maintain optimal tactile sensitivity without compromising protection.

Call for pricing.



Stock Photo - Options May Vary

FEATURES

- Only 0.009" (R-1) thick for decreased finger fatigue
- Rolled and beaded cuff for easy on, easy off
- Thickness: R-1
- Attenuates direct or scattered rays
- 100% tested for pinholes and leaks

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)