

Birkova Altra Gel™ Covered Foam Replacement Spinal Pads: IGC-9015



Stock Photo - Options May Vary

FACTORY NEW

The Birkova Altra Gel™ Covered Foam Spinal Pad Set (IGC-9015) is a replacement pad system designed for Wilson Spine Frames. The pads feature gel molded directly to the foam core and covered with vinyl, eliminating common failure modes such as gel separation and patient skin pinching that occur with standard OEM-style pads.

This pad set provides excellent support and a comfortable surface during spinal procedures. The construction is reusable, disinfectable, and easy to clean. The system is MRI compatible and 100% radiolucent, making it suitable for imaging-compatible surgical environments.

The set includes a 4-piece 1-inch thick Kova Gel™ overlay that attaches via Velcro. The pads are zipper-free and latex-free, offering a cost-effective alternative to standard Wilson Spinal Frame pads with a one-year warranty.

FEATURES

- Cost-effective Wilson Spinal Frames replacement pad set
- MRI compatible & 100% radiolucent
- Gel molded to the foam core
- Covered with vinyl
- Offers great support
- Comfortable surface during procedures
- Zipper-free
- Comes with a 4-pc 1" thick Kova Gel™ overlay set

SPECIFICATIONS

- Replacement pad set for Wilson's Spinal Frame
- Latex-free
- Weight: 35 lbs

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)