

# SchureMed Shoulder Supports | 800-0100

## FACTORY NEW

The SchureMed Shoulder Supports model 800-0100 are surgical table accessories designed to position patients for various surgical procedures on standard operating room tables. These supports feature a broad surface area and clamp securely to accommodate any size patient. Sold in pairs, they provide stable positioning for procedures requiring specific patient angles.

The shoulder supports utilize ball and socket adjustments to enable optimal positioning flexibility for different patient sizes and surgical requirements. They adjust vertically to accommodate X-ray compatible operating table tops, making them versatile for multiple surgical applications.

The supports include deluxe foam pads for patient comfort and are simple to set up and use. Two Schure Socket XL accessories are required for installation.



*Stock Photo - Options May Vary*

## FEATURES

- Simple to set up and easy-to-use
- Comfortable support for Trendelenburg positioning
- Ball and socket adjustments provide optimal positioning for any size patient
- Adjusts vertically to accommodate X-ray tops
- Deluxe foam pads included

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)