

SchureMed Slit Mitt | 807-0054



Stock Photo - Options May Vary

FACTORY NEW

SchureMed Slit Mitt: 807-0054. The Slit Mitt's unique design offers total hand protection while featuring slits for the thumb and fingers, allowing their use during procedures if needed.

This versatile mitt provides excellent protection from radiation, abrasions, and other hazards while maintaining the ability to perform delicate tasks. Its adjustable fit ensures comfort and security, making it ideal for a wide range of surgical and radiographic applications.

The durable material is easy to clean, ensuring long-lasting use in demanding environments.

FEATURES

- Softer lead for increased flexibility
- Reinforced vinyl for strength and ease of cleaning
- Breathable nylon liner to keep hands cooler
- Gloves have finger slits for complete finger movement as well as total hand protection
- Sold as a pair (set of 2)

SPECIFICATIONS

- One size fits all
- 0.5 mm Pb Equiv protection

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)