

SchureMed SpringLoc Clamp | 800-0338



Stock Photo - Options May Vary

FACTORY NEW

The SchureMed SpringLoc Clamp (800-0338) is a side rail clamp designed to quickly and securely attach patient positioning devices and surgical equipment to operating table side rails. The clamp features a flat mounting post and can be positioned anywhere along the rail without requiring table notches, enabling flexible setup for various surgical configurations.

The SpringLoc Clamp utilizes a pull-to-release lever mechanism for quick removal and adjustment during procedures. Its straightforward design allows single-person installation and operation, reducing setup time and improving surgical workflow efficiency.

The clamp provides a strong, secure hold for surgical accessories while maintaining ease of use. Each unit is sold individually for versatile application across different operating room setups and equipment configurations.

FEATURES

- Securely locks accessories with a flat 1" x 1/4" (2.5cm x .6 cm) mounting post
- Attaches easily to the operating table side rail without the need for table notches
- SpringLoc Clamp has a quick & easy pull-to-release lever to remove without the need for a release lever
- Attaches to the surgical table anywhere on the rail
- Single-person installation
- Sold Individually

SPECIFICATIONS

- Length: 3.00" +/- 0.5"
- Width: 2.5" +/- 0.5"
- Depth: 2.25" +/- 0.5"
- Handle Collapsed: 5.5" +/- 0.5"
- Device Weight: 0.5 +/- 0.5 lbs
- The clamp is maintenance-free
- Clamp warranty covers the product from manufacturing defects for two years
- The life of the device is five years under normal use

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)