

SchureMed Surgical Peg Board System | 800-0354



Stock Photo - Options May Vary

FACTORY NEW

The SchureMed Surgical Peg Board System (800-0354) is a complete lateral positioning system designed to optimize surgical positioning and organization. It provides a versatile solution for securing and positioning various surgical accessories during procedures. The system requires two blade-style rail clamps, such as #800-0338, to securely attach the Peg Board to the OR table.

The durable construction ensures reliable performance in demanding surgical environments. The adjustable design allows for customizable positioning to accommodate different surgical needs and procedural requirements.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)