

Steris 3085 SP Surgical Table - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Steris 3085 SP is a versatile surgical table designed for demanding operating room environments. It supports up to 1,000 lbs in standard position and 500 lbs in reverse orientation, making it suitable for a wide range of patients and surgical procedures. The table features full articulation, C-arm compatibility, and durable construction to support efficient surgical workflows.

The 3085 SP includes a wired hand control, battery backup, and hydro-electric brakes for reliable operation. Additional features include return to level functionality, removable and adjustable headrest, reverse orientation patient positioning, and manual kidney lift. North American side rails accommodate most common accessories.

This refurbished unit has been professionally restored and thoroughly tested. It includes a new 2" pad set and wired hand control, providing a complete solution for multi-purpose surgical applications.

FEATURES

- C-Arm Compatible Table Top
- Wired Hand Control
- Battery Operation
- Hydraulic Floor Locks
- Return to Level
- North American Side Rails Fit Most Accessories
- 1,000 lbs Weight Capacity in Normal Position
- Comprehensive Warranty

SPECIFICATIONS

- Weight Capacity: 1000 lbs in Normal Position, 500 lbs in Reverse Orientation
- Table Length (Overall): 76"
- Table Width (Overall): 21 5/8"
- Patient Surface Width: 20"
- Pad Thickness: 2"
- Power Requirements: 120V, 4.5 Amps

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)