

# Wy'East Totalift II Stretcher Chair – Refurbished



Stock Photo – Options May Vary

## FULLY REFURBISHED • CLINICALLY CERTIFIED

The Wy'East Totalift II Stretcher Chair is a versatile patient transport solution that functions as both a chair and stretcher. Its unique dual-purpose design supports patient transport, recovery, and minor procedures across multiple departments and clinical applications.

The unit features pneumatic assist backrest adjustment, hydraulic height control via foot pump, and a sliding patient surface for easy transfers. Central brake and steering systems, fold-under side rails, and an integrated O2 tank well provide comprehensive functionality for clinical environments.

Each unit is completely refurbished to the highest standards and includes a brand-new pad, ensuring comfort and durability. With a 400-pound weight capacity and dependable performance, it delivers a flexible solution for everyday clinical use.

## FEATURES

- Transfer Slide Top for Easy Patient Transfers
- Height Up/Down (Foot Pump at the side of Stretcher)
- Fold-Under Side Rails
- 4 Wheel Locking
- Steering
- 400 lbs Weight Capacity
- Completely Refurbished
- New Pad

## SPECIFICATIONS

- Back-Up/Down (Pneumatic Assist)
- Fold-Under Siderails
- 4 Wheel Locking
- 400 lbs Weight Capacity
- High Position: 36.5"
- Low Position: 21"
- The width of the Pad: is 23"
- Maximum width: 27"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)