

# Steris/Amsco Century 16" Re-Manufactured Sterilizer | V116



Stock Photo - Options May Vary

## FACTORY NEW

Steris/Amsco Century 16" Re-Manufactured Sterilizer model number: V116. Features a Chamber Size of 16" x 16" x 26" and includes a New Reimers Steam Generator, Touch Screen Display, One Rack and Two Shelves. Operating Temperature: 250° to 270°. Cabinet Mounted with a Stainless Steel Enclosure.

Additional Features: The Steri/Amsco V116 has a Foot Pedal Actuated, Automatic Vertical Door, and comes with a Warranty. See for more information.

Note: To view additional models, click the following link, Sterilizers, Disinfectors, & Washers.

## FEATURES

- New Reimers Steam Generator
- One (1) Rack
- Two (2) Shelves
- Touch Screen Display
- Vertical Door: Automatic

## SPECIFICATIONS

- Chamber Size: 16" x 16" x 26"
- Overall Exterior Dimensions (including door lock-wheel): 74.5"H x 30"W x 35.75"D (with door closed)
- Operating Temperature: 250° to 270°
- Cabinet: Stainless Steel
- Steam Generator: Reimers - New

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)