

Steris/Amsco Century 20" Re-manufactured Sterilizer | V120



Stock Photo - Options May Vary

FACTORY NEW

The Steris/Amsco Century 20" Re-manufactured Sterilizer Model Number: V120. The Steris / Amsco Century features a Chamber Size of 20" x 20" x 38", a Touch Screen Display, and includes a New Steam Generator with One Rack and Two Shelves. Cabinet Mounted with a Stainless Steel Enclosure.

Operating Temperature: 250° to 270°. Additional Features: Vertical, Foot Pedal Actuated, Automatic Door and comes with a Warranty.

Note: To view additional models, click the following link, Sterilizers, Disinfectors, & Washers.

FEATURES

- New Reimers Steam Generator
- One (1) Rack
- Two (2) Shelves
- Touch Screen display
- Vertical Door: Automatic

SPECIFICATIONS

- Chamber Size: 20" x 20" x 38"
- Overall Exterior Dimensions (including door lock-wheel): 74.5"H x 32"W x 45.25"D (with door closed)
- Operating Temperature: 250° to 270°
- Cabinet: Stainless Steel
- Steam Generator: Reimers - New

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)