

Stryker 1067 Eye Stretcher – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The Stryker 1067 Eye Stretcher is a specialized surgical stretcher designed for eye, head, and neck procedures. Its tapered head design provides surgeons with optimal access during delicate ophthalmic and cranial operations in surgical and clinical settings.

The stretcher features manual backrest adjustment, hydraulic height control, and Trendelenburg/Reverse Trendelenburg positioning capabilities. It supports patients up to 500 lbs and includes a four-wheel locking system with central brake and steering for reliable mobility and safety.

Additional features include collapsible side rails for patient accessibility, a built-in O2 tank well for respiratory equipment, and a brand-new pad. Each refurbished unit is fully tested and includes a warranty.

FEATURES

- Completely Refurbished
- Manual backrest
- Hydraulic height adjustment
- Trendelenburg/Reverse Trendelenburg
- The tapered head design provides easy access to the surgeon for head, neck, and eye procedures
- 4 Wheel Locking System
- Maximum Weight Capacity of 500 lbs.
- O2 Tank Well

SPECIFICATIONS

- Completely Refurbished
- Manual backrest
- Hydraulic height adjustment
- Trendelenburg/Reverse Trendelenburg
- Collapsible Siderails
- 4 Wheel Locking System
- Maximum Weight Capacity of 500 lbs.
- 6 Months Parts Warranty

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)