

Tuttnauer DS1000 Steam Distiller



Stock Photo - Options May Vary

FACTORY NEW

Tuttnauer DS1000 Steam Distiller. Experience the ease and reliability of the Tuttnauer DS1000 Steam Distiller, your go-to solution for producing high-purity water from any tap source. The DS1000 is simple, efficient, and ideal for clinics, labs, or homes needing clean, contaminant-free water.

The DS1000 uses steam distillation, heating water in a stainless-steel boiler to 212°F (100°C), the boiling point. The high temperature kills bacteria, viruses, and other microorganisms, ensuring the water is as clean and pure as possible.

Furthermore, once the water vapor condenses, it collects in a clean container, leaving impurities behind. Operating the DS1000 couldn't be easier.

FEATURES

- Produces 1 gallon of 100% steam distilled water in 4 hours (Up to 6 gallons per day)
- Reduces 99% of total dissolved solids (TDS)
- Dripleless stop-and-serve feature
- Automatic shutoff after each 1 Gallon cycle
- Effectively removes arsenic, nitrates, sodium, chlorine, bacteria and viruses

SPECIFICATIONS

- 800 watts of power, 120 Volt 60 Hz, fits standard outlet
- Item Weight 11.43 pounds
- Dimensions: 15.25" H x 15" W x 9.75" D

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)