

# Tuttnauer 1730M Valueklave Manual Autoclave



Stock Photo - Options May Vary

## FACTORY NEW

The Tuttnauer 1730M Valueklave Manual Autoclave is a compact sterilization unit designed for small offices or as a secondary autoclave for larger practices. It features simple controls and an easy-to-close door for straightforward operation. The unit delivers fast cycle times and heated drying for quick turnaround of wrapped or unwrapped instruments.

The 1730M incorporates an electropolished chamber and door for long-term durability, with a front-mounted drain valve for convenient water reservoir draining. Dual safety mechanisms protect against overheating, and a double safety locking device prevents door opening while the chamber is pressurized. The autoclave includes international certifications and comes with a 1-year parts and labor warranty.

## FEATURES

- Fast Cycle Time, 11 min from a hot start\* (Standard Unwrapped Cycle Time)
- Automatic Shut Off at the End of Both the Sterilization and Dry Cycles
- A Long Life Electro Polished Chamber and Door
- Double Safety Locking Device Prevents Door from Opening while the Chamber is Pressurized
- The Drain Valve is Located on the Front, Allowing for Quick and Easy Draining of Water Reservoir
- Dual Safety to Protect Against Overheating
- International Certifications
- 1 Year Parts and Labor Warranty

## SPECIFICATIONS

- Chamber Dimensions: 7" x 13"
- Chamber Volume: 2 Gallons
- Overall Dimensions: 17.9"D x 17.4"W x 12"
- Standard Cassettes Capacity: 2 Half
- Number of Trays: 3
- Standard Unwrapped Cycle Time: Cold: 16 Min, Hot: 11 Min
- Volts: 120V
- Power: 1350W

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)