

Tuttnauer 2340M Manual Autoclave



Stock Photo - Options May Vary

FACTORY NEW

The Tuttnauer 2340M Manual Autoclave. The Tuttnauer 2340M Manual Autoclave is a trusted, high-performance sterilizer designed to meet the rigorous demands of medical, dental, and laboratory environments. The 2340M's durable design and bright, user-friendly control panel make every sterilization cycle simple and efficient.

Engineered for safety and efficiency, this manual autoclave includes an automatic shut-off at the end of both the sterilization and drying cycles, helping reduce energy use and staff oversight. A double safety locking mechanism ensures the chamber cannot be opened while pressurized, providing critical protection for operators.

The compact and reliable Tuttnauer 2340M fits easily into clinics, labs, and centers.

FEATURES

- Automatic Shut Off at the end of Both the Sterilization and Dry Cycles
- A Long Life Electro Polished Chamber and Door
- Double Safety Locking Device Prevents Door from Opening while the Chamber is Pressurized
- The Drain Valve is Located on the Front, allowing for Quick and Easy Draining of Water Reservoir
- Dual Safety Thermostat to Protect Against Overheating

SPECIFICATIONS

- Chamber Dimensions: 9" x 18"
- Chamber Volume: 5 Gallons
- Overall Dimensions: 21.5"D x 20"H x 14.4"H
- Standard Cassettes Capacity: 2 Full and 2 Half
- Tray dimensions: 16.3"D x 6.7"W x 0.8"H

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)