

Tuttnauer 2540M Manual Autoclave



Stock Photo - Options May Vary

FACTORY NEW

The Tuttnauer 2540M Manual Autoclave. Tuttnauer manual autoclaves set the standard for steam sterilization. Our basic manual autoclave will satisfy all of your sterilization needs without compromising on quality, safety, or reliability.

The 2540M complements any healthcare facility and includes 4 full-size trays. Additional 2540M Features: The design team created the bright, easy-to-use panel with the operator in mind. Additionally, the exceptionally deep chamber accommodates 6 standard cassettes—3 full and 3 half.

To ensure safety, the system features automatic shut-off at the end of both the Sterilization and Dry Cycles. Furthermore, a Double Safety Locking Device prevents the door from opening while the chamber is pressurized.

FEATURES

- Automatic Shut Off at the end of Both the Sterilization and Dry Cycles
- Long Life Electro Polished Chamber and Door
- Double Safety Locking Device Prevents Door from Opening while the Chamber is Pressurized
- Drain Valve is Located on the Front, Allowing for Quick and Easy Draining of the Water Reservoir
- Dual Safety Thermostat to Protect Against Overheating

SPECIFICATIONS

- Chamber Dimensions: 10" x 19"
- Chamber Volume: 6 Gallons
- Overall Dimensions: 21.5"D x 20"W x 14.4"H
- Standard Cassettes Capacity: 3 Full and 3 Half
- Tray Dimensions: 16.3"D x 6.7"W x 0.8"H

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)