


Tuttnauer 2540MK Manual Kwiklave



Stock Photo - Options May Vary

FACTORY NEW

The Tuttnauer 2540MK Manual Kwiklave is available through Venture Medical. This sterilizer offers all the proven, dependable features of our manual autoclave, with the added benefit of reducing the total heat-up time.

This autoclave is the perfect choice for a quick instrument turnaround or a high-volume office. Additional 2540MK Features: Speed and reliability are accomplished without sacrificing load size. Automatic Shut Off at the End of Both the Sterilization and Dry Cycles, along with a Long Life, Electro-Polished Chamber and Door. Double Safety Locking Device and Dual Safety Thermostat to Protect against Overheating. Additionally, this unit has International Certifications and comes with a 1-Year Parts and Labor Warranty.

FEATURES

- Automatic Shut Off at the End of Both the Sterilization and Dry Cycles
- A Long Life Electro Polished Chamber and Door
- Double Safety Locking Device Prevents Door from Opening while the Chamber is Pressurized
- Drain Valve is Located on the Front, Allowing for Quick and Easy Draining of the Water Reservoir
- Dual Safety Thermostat to Protect Against Overheating

SPECIFICATIONS

- Chamber Dimensions: 10" x 19"
- Chamber Volume: 6 Gallons
- Overall Dimensions: 21.5"D x 20"W x 14.4"H
- Standard Cassettes Capacity: 3 Full and 3 Half
- Tray Dimensions: 16.3"D x 6.7"W x 0.8"H

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)