

Tuttnauer 3870M Manual Autoclave



Stock Photo - Options May Vary

FACTORY NEW

Tuttnauer 3870M Manual Autoclave. The Tuttnauer 3870M Manual Autoclave is the largest tabletop autoclave in Tuttnauer's manual lineup, making it the ideal sterilization solution for multi-doctor practices, busy clinics, and surgery centers. With a large 15" x 30" chamber and 85-liter capacity, this autoclave fits 15 cassettes or long instruments easily.

Built for power, precision, and safety, the 3870M features automatic shut-off at the end of both sterilization and drying cycles, ensuring user convenience and energy efficiency. The double-locking system and dual thermostats prevent door opening under pressure and protect against overheating for safer operation.

Certified to meet international standards, the Tuttnauer 3870M delivers dependable performance with every cycle.

FEATURES

- Chamber Volume: 85 Liters
- Standard Cassette Capacity: 15 Full
- Number of Trays: 2
- Automatic shut-off at the end of both the sterilization and dry cycles
- A long-life electro-polished chamber and door

SPECIFICATIONS

- Chamber Dimensions: 15" x 30"
- Chamber Volume: 22 Gallons
- Overall Dimensions: 34.5"D x 26"W x 20.6"H
- Standard Cassettes Capacity: 15 Full
- Tray Dimensions: 26.6"D x 11.3"W x 1"H

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)