

Tuttnauer 7000 Series Distillers



Stock Photo - Options May Vary

FACTORY NEW

The Tuttnauer 7000 Series is a high-capacity countertop water distiller designed for medical, dental, and laboratory facilities. It produces 99% TDS-free water, offering a reliable, cost-effective alternative to bottled water for autoclaves, lab instruments, and sensitive procedures.

The unit features direct water line connection, enabling continuous, hands-free operation that minimizes downtime and ensures a steady supply of distilled water. The 3.5-gallon reservoir supports front or side spout mounting for flexible installation in various facility layouts.

The distiller is designed for ease of use, cleaning, and maintenance, making it ideal for busy environments that require pure water on a daily basis.

FEATURES

- Directly connected to the water supply
- Allows for continuous operation
- Front or side mounting of the spout
- 3.5 Gallon System

SPECIFICATIONS

- Model Type: Countertop
- Weight: 30 lbs
- Width: 11"
- Depth: 17"
- Height: 17"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)