

Tuttnauer 7000-12 Series Distiller



Stock Photo - Options May Vary

FACTORY NEW

The Tuttnauer 7000-12 Large Capacity Distiller delivers a reliable, high-volume solution for facilities requiring a steady supply of pure water. This distiller produces water 99% free of total dissolved solids (TDS) and eliminates the cost, waste, and hassle of bottled water.

With a 12-gallon capacity, the 7000-12 is built for continuous performance. The unit connects to your water line for automatic refilling and nonstop use in busy environments.

Designed as a floor-standing model, the Tuttnauer 7000-12 maximizes output while minimizing operator effort. Its efficient design and durable construction make it a trusted solution for medical, dental, surgical, and laboratory settings.

FEATURES

- Floor Model
- Directly connected to the water supply
- Allows for continuous operation
- Front or side mounting of the spout
- 12 Gallon System

SPECIFICATIONS

- Model Type: Floor Model
- Weight: 44 lbs
- Width: 11"
- Depth: 17"
- Height: 42"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)