

Tuttnauer 7000-8 Series Distiller



Stock Photo - Options May Vary

FACTORY NEW

The Tuttnauer 7000-8 Large Capacity Distiller produces 99% TDS-free water, reducing bottled water costs while reliably supplying clean water for autoclaves and laboratory equipment. With an 8-gallon capacity, the 7000-8 is ideal for medium- to high-volume facilities in medical, dental, and research settings.

The distiller connects directly to the water line and operates automatically, enabling continuous operation without manual intervention. Its floor-standing design integrates practical space efficiency into busy clinical environments.

The reservoir accommodates front or side mounting of the spout, offering flexible installation options based on facility layout. The unit is straightforward to use, clean, and maintain, providing long-term value for facilities requiring purified water.

FEATURES

- Directly connected to the water supply
- Allows for continuous operation
- Front or side mounting of the spout
- 8 Gallon System

SPECIFICATIONS

- Model Type: Floor Model
- Weight: 37 lbs
- Width: 11"
- Depth: 17"
- Height: 36"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)