

Tuttnauer Autoclave 3870EAP Fully Automatic



Stock Photo - Options May Vary

FACTORY NEW

Tuttnauer Autoclave 3870EAP Fully Automatic. The Tuttnauer Autoclave 3870EA Fully Automatic is our largest tabletop autoclave, ideal for high-volume sterilization. It features a closed-door active drying system, which keeps instruments sterile and ensures efficient drying for packs and pouches. The Tuttnauer 3870EAP also includes an internal printer.

This autoclave is perfect for multi-doctor practices, surgery centers, and busy clinics needing fast, reliable sterilization. The 15" x 30" chamber holds longer instruments or large quantities of smaller tools. You can fit up to 15 full-size IMS cassettes at once.

Additionally, it's a fully automatic system that streamlines sterilization from start to finish, saving valuable staff time.

FEATURES

- Closed-door active drying system to maintain sterility and ensure efficient drying of packs and pouches
- Hepa filter 0.2 µm air filter (replaceable)
- Rapid cycles optimize the availability of ready-to-use sterile instruments
- Easy-to-use, single touch control system increases productivity
- 3 Pre-set sterilization programs

SPECIFICATIONS

- Model Number: 3870
- Chamber Dimensions; 15" x 30"
- Chamber Volume: 22 Gallons
- Standard Cassettes Capacity: 15 Full
- Number of Trays: 2

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)