

Tuttnauer CSU1 Ultrasonic Cleaner



Stock Photo - Options May Vary

FACTORY NEW

Tuttnauer CSU1 Ultrasonic Cleaner - Manufactured with only the highest quality components. These machines offer 60 kHz of gentle, effective cleaning power, safely removing surface contamination without damaging instruments. Furthermore, contemporary design and quiet operation help the CSU1 complement any office environment.

Clean & Simple Ultrasonic Systems come in two sizes with accessories to match your specific cleaning needs. Additionally, the 1-gallon system (shown) holds a standard selection of instruments or 2 half-size cassettes. Do you have questions?

To view a wide variety of Autoclaves & Equipment, click the following link to browse automatic autoclaves and distillers to manual autoclaves, and ultrasonic cleaners such as the CSU1.

FEATURES

- High efficiency 60 kHz cavitation for a thorough and quieter cleaning
- Durable vinyl clad steel housing which is easy to clean and resists scratching
- Tank Capacity: 1 Gallon
- Stainless Steel Tank with Cover

SPECIFICATIONS

- Tank Capacity: 1 Gallon
- Internal Tank Depth: 5-3/8"
- Internal Tank Width: 9-3/8"
- Internal Tank Height: 4"
- External Depth: 9-7/8"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)