

Tuttnauer CSU3 Ultrasonic Cleaner



Stock Photo - Options May Vary

FACTORY NEW

Tuttnauer CSU3 Ultrasonic Cleaner is manufactured with only the highest quality components. These machines offer 60 kHz quiet cleaning, gently removing surface contamination without damaging instruments. Furthermore, contemporary design and quiet operation complement any office environment.

CSU3 units come in two sizes with a wide range of accessories to fit your specific ultrasonic cleaning needs. Furthermore, the 3-gallon system (shown) holds a standard selection of instruments or full-size cassettes. Do you have questions?

To view a wide variety of Autoclaves & Equipment, click the following link to browse automatic autoclaves and distillers to manual autoclaves, and ultrasonic cleaners like the CSU3. Venture Medical's Complete Catalog of Autoclaves.

FEATURES

- High efficiency 60 kHz cavitation for a thorough and quieter cleaning
- Durable vinyl clad steel housing is easy to clean and resists scratching
- 3 Gallon
- Stainless Steel Tank with Cover

SPECIFICATIONS

- Internal Tank Width: 15-1/2"
- Internal Tank Height: 6"
- Internal Tank Depth: 9"
- External Width: 19"
- External Height: 13-1/2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)