

Tuttnauer CSU3H Ultrasonic Cleaner with Heater



Stock Photo - Options May Vary

FACTORY NEW

Tuttnauer CSU3H Ultrasonic Cleaner with Heater. Clean & Simple Ultrasonic Cleaning Systems with Heaters are manufactured with only the highest quality components. These machines deliver 60 kHz of quiet cleaning power, gently cleaning instruments while effectively removing all surface contamination.

Additional CSU3H Features: The contemporary design and quiet operation complement any office environment. Furthermore, Clean & Simple Ultrasonic Systems are available in 2 sizes with a full selection of accessories to meet your specific needs. Additionally, the 3-gallon system (shown) holds a standard selection of instruments or full-size cassettes.

Tuttnauer also offers the smaller 1-gallon system, which holds 2 half-size cassettes. Do you have questions?

FEATURES

- High efficiency 60 kHz cavitation for a thorough and quieter cleaning
- Heater
- Durable vinyl clad steel housing is easy to clean and resists scratching
- 3 Gallon
- Stainless Steel Tank with Cover

SPECIFICATIONS

- Internal Tank Width: 15-1/2"
- Internal Tank Height: 6"
- Internal Tank Depth: 9"
- External Width: 19"
- External Height: 13-1/2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)