

Tuttnauer EZ10 Fully Automatic Autoclave



Stock Photo - Options May Vary

FACTORY NEW

Tuttnauer EZ10 Fully Automatic Autoclave. The Tuttnauer EZ10 makes sterilization easy with one-touch operation—perfect for busy medical and dental offices. Additionally, it automatically fills, sterilizes, exhausts, and dries with the push of a button.

This autoclave saves time, reduces effort, and delivers reliable results every cycle. Additionally, the Tuttnauer EZ10 Autoclave features a user-friendly control panel and intuitive design for effortless operation.

Key features include one-touch operation, automatic shut-off, a double safety locking system, and a dual safety thermostat for enhanced protection. The EZ10 is internationally certified and comes with a 2-year parts and labor warranty, offering long-term value and peace of mind.

FEATURES

- Closed-door active drying system to maintain sterility and ensure efficient drying of packs and pouches
- A High-efficiency air pump circulates hot air after the sterilization cycle for fast drying
- Hepa filter 0.2 µm air filter (replaceable)
- Pre-loaded cycles for unwrapped instruments, wrapped/packs, glassware, and additional drying cycle
- All program parameters can be customized and stored

SPECIFICATIONS

- Chamber Dimensions: 10" x 19"
- Chamber Volume: 6 Gallons
- Overall Dimensions: 21.5"D x 20"W x 14.4" H
- Standard Cassettes Capacity: 3 full and 3 half
- Tray dimensions: 16.3"D x 6.7"W x 0.8"H

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)