

# Tuttnauer EZ10K Fully Automatic Kwiklave



Stock Photo - Options May Vary

## FACTORY NEW

Tuttnauer EZ10K Fully Automatic Kwiklave. The Tuttnauer EZ10K delivers fast, reliable sterilization for high-demand practices. It reduces heat-up time significantly, making it ideal for quick instrument turnaround or busy, high-volume offices. Built on the trusted EZ series platform, it combines speed, safety, and performance in one compact unit.

This autoclave features a powerful, single-touch control system that streamlines operation. Staff can easily choose between three pre-set sterilization programs, increasing workflow efficiency. The EZ10K also allows full cycle customization, offering flexibility for various instruments and load types.

Despite its speed, the EZ10K does not compromise on capacity.

## FEATURES

- Reduces Cleaning Cycle Time by 50%
- Closed-door active drying system to maintain sterility and ensure efficient drying of packs and pouches
- Hepa filter 0.2 µm air filter (replaceable)
- Rapid cycles optimize the availability of ready-to-use sterile instruments
- Easy-to-use, single touch control system increases productivity

## SPECIFICATIONS

- Model Number: EZ10K
- Chamber Volume: 6 Gallons
- Chamber Dimensions: 10" x 19"
- Standard Cassettes Capacity: 3 full and 3 half
- Number of Trays: 4

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)