

Tuttnauer EZ9 Fully Automatic Autoclave



Stock Photo - Options May Vary

FACTORY NEW

The Tuttnauer EZ9 is a fully automatic autoclave designed for clinics, dental offices, and medical practices requiring reliable sterilization. One-touch operation automatically fills, sterilizes, exhausts, and dries with a single button press. The closed-door active drying system maintains sterility while efficiently drying both packs and pouches using a high-efficiency air pump that circulates hot air after sterilization.

The EZ9 features a HEPA 0.2 μm replaceable air filter that filters air entering the chamber during the drying process, enhancing safety and cleanliness. Users can set and save custom cycle settings to fit specific sterilization needs, and the autoclave automatically shuts off after sterilizing and drying to conserve energy.

The unit includes audible and visual alerts if a cycle is interrupted, enabling staff to respond quickly. The EZ9 holds international certifications and comes with a 2-year parts and labor warranty.

FEATURES

- Closed-door active drying system to maintain sterility and ensure efficient drying of packs and pouches
- Hepa filter 0.2 μm air filter (replaceable)
- Pre-loaded cycles for unwrapped instruments, wrapped/packs, glassware, and additional drying cycle
- All program parameters can be customized and stored
- Audible and visual cycle interruption alert

SPECIFICATIONS

- Chamber Dimensions: 9" x 18"
- Chamber Volume: 5 Gallons
- Standard Cassettes Capacity: 2 Full and 2 Half
- Voltage: 120V
- Power: 1400W

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)