

Tuttnauer EZ9Plus Fully Automatic Autoclave



Stock Photo - Options May Vary

FACTORY NEW

Tuttnauer EZ9Plus - Fully Automatic Autoclave. The EZ9Plus meets ANSI/AAMI ST55 sterilization standards, helping practitioners manage challenging workloads and regulatory requirements efficiently. The newly designed, ergonomically friendly EZGlide door features a sleek, modern look and effortless opening and closing.

The EZ9Plus Fully Automatic Autoclave offers an easy-to-read EZView Multi-Color Display that shows real-time cycle updates clearly. It displays key details like temperature, pressure, and cycle steps, with each stage shown using a color wheel and simple text.

Plus, the closed-door drying feature keeps items sterile while drying, making the process safer and more efficient.

FEATURES

- Newly designed, ergonomically friendly EZGlide door features a sleek modern look and effortless opening and closing
- Innovative EZView Multi-Color Display offers real-time cycle information that is easy to read, even from a distance. The EZView displays cycle parameters, temperature, and pressure in the chamber. During a cycle, each stage is displayed with an active color wheel and text
- Closed-door drying, meaning the door remains closed through the dry cycle to maintain sterility and ensure efficient drying. A high-efficiency air pump circulates filtered air through the chamber for fast drying.
- The 0.2 μ HEPA filter is replaceable
- Comprehensive 2-year parts & labor warranty and 10-year chamber warranty

SPECIFICATIONS

- Chamber Dimensions: 9" x 19.8
- Volume: 5.2 Gallons
- Overall Dimensions: 24.8" D x 19.1" W x 15.1" H
- Standard Cassette Capacity: 2 Full & 2 Half (loaded vertically)
- Tray Dimensions: 16.3" D x 6.7" x 0.6" H

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)