

UMF 24" Base Cabinet | 6014



Stock Photo - Options May Vary

FACTORY NEW

UMF 24" Base Cabinet: 6014. UMF modular casework allows you to be creative in setting up an exam room. Begin the process by selecting from the many different base cabinet configurations offered. All cabinets are made of All-Steel construction.

All cabinets can be fitted with optional locks. Custom sizes and configurations are available. The UMF 24" Base Cabinet Features: 20 gauge steel construction and a high-pressure laminate countertop with a 4" (10.2 cm) backsplash. Full-extension roller bearing drawer glides along with Commercial grade, adjustable; spring assisted door hinges.

The drawer and door fronts are interchangeable. The 6014 is available in all UMF Medical colors and has a 3-year limited warranty.

FEATURES

- 20 gauge steel construction
- High-pressure laminate countertop with 4" backsplash
- Full extension roller bearing drawer glides
- Commercial grade, adjustable, spring-assisted door hinges
- Interchangeable drawer and door fronts

SPECIFICATIONS

- 24" wide
- 34.5" high
- 18" deep
- 24" wide
- 35.875" high - does not include the height of the backsplash

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)