

UMF 24" Wall Cabinet | 6074



Stock Photo - Options May Vary

FACTORY NEW

UMF 24" Wall Cabinet: 6074. UMF commits to providing excellence to healthcare providers, helping them serve their patients effectively. For nearly 60 years, our customers have relied upon the excellent construction, durability, and reliability of UMF Medical products. Furthermore, each product promotes excellence, efficiency, and service in delivering healthcare services.

UMF 24" Wall Cabinet 6074 Features: UMF designs wall cabinets to complement the modular base cabinets. They are 20" high and 11" deep. The cabinets feature a welded "J" hook and include wall hanger strips with holes at 1/2" intervals for easy installation and leveling.

All wall cabinets feature 20 gauge steel construction and have commercial-grade, easily adjustable, spring-assisted, concealed door hinges.

FEATURES

- 20 gauge steel construction
- Commercial grade, adjustable, spring-assisted door hinges
- Adjustable shelf
- Interchangeable door fronts
- Wall mounting bracket

SPECIFICATIONS

- 24" Wide
- 20" High
- 11" Deep

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)