

UMF Bariatric Power Treatment Table | 5005



Stock Photo - Options May Vary

FACTORY NEW

UMF 5005 Bariatric Power Treatment Table. UMF is committed to providing excellence for healthcare providers in service to their patients. For nearly 60 years, our customers have relied upon the excellent construction, durability, and reliability of UMF Medical products.

Moreover, we build each product to promote excellence, efficiency, and exceptional healthcare service delivery. The UMF Medical 5005 Bariatric Treatment Table provides unmatched capacity and comfort in a stable, heavy-duty all-steel design. In addition, the 5005 has an 800-lb.

(362.9kg) patient weight capacity and a low 19" seat height. It also features a fully articulating power back section as well as power hi-lo control.

FEATURES

- Ultra comfort seamless top with antimicrobial guard and built-in pillow
- ADA Compliant with 19" (48.3cm) low access
- All-steel heavy-duty construction
- Power Hi-Lo
- Two-function hand control

SPECIFICATIONS

- Height: Min. 19" (48.3cm) at seat Max. 35" (88.9cm) at the seat
- Width: 30" (76.2cm)
- Length: 70" (177.8cm) to 80" (203.2cm) with leg section extended
- Storage drawer: 16" D x 19" W x 4.5" H (40.6cm D x 48.3cm W x 11.4cm H)
- Electrical: 120 VAC, 60 HZ or 220/240 VAC, 60HZ

**Brochure and images are a general overview of the product. Actual product options may vary.*