

UMF Medical Power ENT/Phlebotomy Chair | 8678

FACTORY NEW

UMF 8678 Medical Power ENT/Phlebotomy Chair. UMF is committed to providing excellence for healthcare providers in service to their patients. Furthermore, for nearly 60 years, our customers have relied upon the excellent construction, durability, and reliability of UMF Medical products.

Additionally, each product is built to promote excellence, efficiency, and service in the delivery of healthcare services. The UMF 8678 Power ENT/Phlebotomy Chair offers comfort and an electro-mechanical actuator. It features fully articulating arms for easy access.

The chair's 19" low access is ADA-compliant, aiding patients with limited mobility. The two-function hand control adjusts seat height and backrest angle.



Stock Photo - Options May Vary

FEATURES

- ADA-compliant, 19" low access height
- 24" low access height
- 375-pound patient weight capacity
- Power Hi-Lo
- Power backrest

SPECIFICATIONS

- Height: Min. 19" (48.3 cm) at seat
- Max. 35" (88.9 cm) at seat
- Overall Width: 31.5" (80.0 cm)
- Seat Width: 18" (45.7 cm)
- Length (fully extended): 67" (170.2 cm)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)