

UMF Power Exam Chair | 4011



Stock Photo - Options May Vary

FACTORY NEW

UMF 4011 Power Exam Chair. Ultra Procedure Chair with One-Touch Patient Positioning® System Hand Control For Procedures. The Power4011 is a product of extensive research and clinical feedback.

Furthermore, it was developed to improve patient care and streamline healthcare providers' work. Uses extend to procedures performed by OB/GYN and Oncology professionals. Designed for procedures, the power4011 meets ADA and U.S.

Access Board regulations and accommodates the needs of most patients. UMF manufactures all medical products in the U.S.A.

FEATURES

- 650 lb. (294.84kg) patient weight capacity
- Low access ADA height 18.5" (±0.5") 46.9cm (±1.27cm)
- 4 power functions: hi-lo, back, tilt, leg lift & extension
- Hand control with one-touch patient positioning® system, 2 programmable positions & return-to-home
- Removable, seamless upholstery with Prefixx® Protective Finish

SPECIFICATIONS

- Height: min. 18.5" (± .5") 46.9cm (±1.27cm) at seat max. 37" at seat (93.9cm)
- Seat width: 28" (71.12cm)
- Back width: 26.5" (67.31cm)
- Length: 73"-80" with foot section extended (185.42cm - 203.2cm)
- Leg rest: 28" W x 11.5" D (71.12cm W x 29.21cm D)

**Brochure and images are a general overview of the product. Actual product options may vary.*