

UMF Power Podiatry Table w/Accessories | 5016



Stock Photo - Options May Vary

FACTORY NEW

UMF 5016 Power Podiatry Table w/Accessories. UMF is committed to providing excellence for healthcare providers in service to their patients. For nearly 60 years, our customers have relied upon the excellent construction, durability, and reliability of UMF Medical products. UMF builds each product to promote excellence, efficiency, and outstanding healthcare service delivery.

UMF 5016 Power Podiatry Table w/Accessories Features: The ADA-Compliant UMF Medical 5016 Power Podiatry Chair redefines easy patient access with its fully articulating arms, adjustable leg section, and 19" (48 cm) low access seat height. With its three-function memory position foot control, you'll have full control of seat height, backrest angle, and seat tilt angle.

FEATURES

- ADA Compliant - 19" (48.3 cm) low access seat height
- UMF Medical Ultra Comfort Seamless Top
- ADA Compliant
- 19" (48 cm) low-access seat height
- Power Hi-Lo

SPECIFICATIONS

- Height: Min. 19" (48.3 cm) at seat
- Max. 35" (88.9cm) at seat
- Width: 26.5" (67.3 cm)
- Length upright: 50" to 59" (127.0 cm to 149.9 cm)
- Length flat: 74" to 83" (188.0 cm to 210.8 cm)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)