

UMF Signature Series Exam Table | 5240-145



Stock Photo - Options May Vary

FACTORY NEW

UMF Signature Series Exam Table: 5240-145. UMF is committed to providing excellence for healthcare providers in service to their patients. For nearly 60 years, our customers have relied upon the excellent construction, durability, and reliability of UMF Medical products.

Each product promotes excellence, efficiency, and service in the delivery of healthcare services. The UMF Medical 5240-145 Signature Series Exam Table Features: The UMF Medical 5240 Signature Series exam table maintains the tradition of strength and quality of the original UMF Medical Signature Series but with an updated look.

FEATURES

- Removable UMF Medical Ultra Comfort Seamless Top
- Articulating back section with adjustable backrest
- One-touch, Easy-glide comfort stirrups
- Slide out leg section with removable upholstered leg rest
- Removable plastic drip pan 9" x 9" / 22.9cm x 22.9cm)

SPECIFICATIONS

- Height: 31.8"
- Length: 58.5" - 77.5"
- Width: 27"
- Stirrups: 3 lateral positions, infinitely adjustable over the length of the bar

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)