

UMF Signature Series Exam Table | 5250



Stock Photo - Options May Vary

FACTORY NEW

UMF Signature Series Exam Table: Model Number 5250. UMF commits to providing excellence for healthcare providers to serve their patients.

For nearly 60 years, our customers have relied upon the excellent construction, durability, and reliability of UMF Medical products. Furthermore, each product is built to promote excellence, efficiency, and service in the delivery of healthcare services.

The UMF 5250 Signature Series Exam Table Features: The 5250 continues the UMF Medical Signature Series tradition of strength and quality, with a fresh, updated look. Reversible, integrated molded drawers provide flexibility in room layout and eliminate the need for choosing a right or left-handed table.

FEATURES

- Removable UMF Medical Ultra Comfort Seamless Top
- Articulating back section with gas spring assist
- One-touch, Easy-glide comfort stirrups
- Slide out leg section with removable upholstered leg rest
- Removable plastic drip pan (9" x 9" / 22.9 cm x 22.9 cm)

SPECIFICATIONS

- Height: 31.8"
- Length: 58.5" - 77.5"
- Width: 27"
- Stirrups: 3 lateral positions, infinitely adjustable over the length of the bar

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)