

ValleyLab Force 1C ESU – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The ValleyLab Force 1C ESU is a versatile, compact unit that is great for Endo. Features Include Enhanced Safety with REM Return Electrode Monitoring, which helps Prevent Electrode Burns. The max cutting power is 200 watts.

The unit includes a monopolar footswitch and a user manual. An optional cart and bipolar footswitch are available. Are you looking for a specific Electro-surgical Unit or Accessory? Click the following link to view our complete catalog of New or Completely Refurbished Electro-Surgical Units. Have Questions?

Contact us today for questions or available options to fit your specific needs!

FEATURES

- Versatile, Compact Unit, Great for Endo
- 200W Cut / 75W Coag/ 50W Bipolar
- Enhanced Safety with REM Return Electrode Monitoring
- Audible Alarms
- 3 Preset Blend Modes

SPECIFICATIONS

- 200W Cut
- 75W Coag
- 50W Bipolar

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com