

# Champion Phlebotomy Chair with Drawer | The Spirit 2570



Stock Photo - Options May Vary

## FACTORY NEW

The Champion Phlebotomy Chair Spirit 2570 is a high-quality blood drawing chair designed for use in clinics, labs, outpatient centers, and hospitals. It features two infinitely adjustable armrests, a durable one-piece resin seat and backrest, and non-marring leveling feet for stability on various floor surfaces. The flip-up contoured urethane armrest provides comfort and easy access while maintaining safety for patients during phlebotomy procedures.

The Spirit 2570 is built with a high-strength steel frame and offers adjustable arm height to accommodate different patient needs and operator preferences. The chair includes convenient drawer storage for supplies and equipment. With a 300-pound weight capacity, it is suitable for a wide range of patients.

Champion has established a reputation for delivering exceptional quality and functionality in medical seating. The Spirit 2570 is engineered to provide years of reliable, trouble-free service in any healthcare setting.

## FEATURES

- High-Strength Steel Frame
- Infinitely Adjustable Height
- Flip-Up Armrest
- Convenient Drawer Storage Cabinet

## SPECIFICATIONS

- Seat Height: 20.5"
- Seat Width: 18"
- Seat Depth: 15"
- Overall Height: 36"
- Overall Width: 35.25"
- Armrest Adjustable Height: 26.5"-32.5"
- Weight Capacity: 300 lbs
- Model Number: 2570

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)