

Winco Blood Drawing Chair with Padded Vinyl | 2573



Stock Photo - Options May Vary

FACTORY NEW

The Winco Blood Drawing Chair model 2573 is designed for patient comfort during phlebotomy procedures. The chair features a padded seat and back constructed with 1.25" medium-density foam, supported by a high-strength steel frame with non-slip rubber feet for stability and safety.

The chair includes comfortable padded armrests available in either a standard flip-up urethane or optional padded vinyl pivot configuration. A neutral gray powder coat finish provides durability and a professional appearance suitable for any clinical setting.

The chair is available in 10 different colors to match facility décor preferences. Its robust construction and thoughtful design make it an ideal choice for blood draw stations and phlebotomy centers.

FEATURES

- Padded seat and back
- High-strength steel frame
- Non-slip rubber feet
- Comes in 10 different colors to match any decor
- A neutral gray powder coat provides a durable finish
- Cushioned with 1.25" medium-density foam
- Comfortable padded armrests

SPECIFICATIONS

- Seat Height: 22.5"
- Seat Width: 18.75"
- Overall Height: 44.5"
- Overall Width: 29"
- Armrest Adj. Height: 28.5"-38.5"
- Overall Depth: 27"
- Weight: 48 lbs
- Weight Capacity: 300 lbs

**Brochure and images are a general overview of the product. Actual product options may vary.*