

Winco Blood Drawing Chair | 2571

FACTORY NEW

Winco Blood Drawing Chair: 2571. Value and quality at an affordable price. Economical, yet designed with durability and longevity in mind. The high-quality features of Winco's Blood Drawing Chairs will provide you with years of service.



Stock Photo - Options May Vary

Winco's 2571 features infinitely adjustable armrests, a durable plastic resin seat, and non-marring leveling feet. Additional 2571 Features: The flip-up contoured urethane armrest provides comfort, easy access, and added safety for the user. If a cabinet is required, please see our 2570 Blood Drawing Chair.

Neutral gray powder coating provides a durable finish. Please Note: Winco's 2570 & 2571 Blood Drawing Chair (plastic seat and back) comes in Black only.

FEATURES

- Durable plastic Seat/Back
- Sturdy Leveling Feet
- High-strength steel frame
- Infinitely adjustable height
- Flip-up armrest

SPECIFICATIONS

- Seat Height: 20.5"
- Seat Width: 18"
- Seat Depth: 15"
- Overall Height: 36"
- Overall Width: 29"

**Brochure and images are a general overview of the product. Actual product options may vary.*