

# Champion T3 Procedure Chair



Stock Photo - Options May Vary

## FACTORY NEW

The Champion T3 Procedure Chair provides a versatile diagnostic seating solution for modified barium swallow studies and video fluoroscopy procedures. The seat rotates 90 degrees to capture clear lateral and A/P imaging without repositioning the patient, improving safety and streamlining clinical workflows. A radiolucent back ensures high-quality imaging results during every diagnostic session.

Powered back, leg, and height adjustments allow staff to transition from a procedure chair to a transport device using simple push-button pendant control. Swing and tuck arm rails create extra space for tight radio frequency environments. The compact wheelbase navigates easily through narrow hallways and standard elevators, while the heavy-duty frame supports up to 350 pounds.

A removable lithium-ion battery powers the unit for extended use in busy clinics. The fold-down push bar provides caregivers with a secure grip for smooth transport throughout the hospital. The T3 is available in three models: T3-B, T3-FB, and T3-SB.

## FEATURES

- 90° Locking Rotation: Rotate the seat to capture lateral and A/P views without repositioning the patient.
- Radiolucent Backrest: High-clarity materials ensure superior imaging quality for barium swallow and fluoroscopy studies.
- One Patient, One Surface: Streamline workflows by using one device for transport, procedures, and recovery.
- Hand Control for Powered Articulation: Adjust the height, back, and leg sections via a simple push-button hand control pendant.
- Swing-Away Side Rails: Tuck the armrests away to provide maximum clearance for imaging equipment.
- Lithium-Ion Battery: The removable battery system allows for cordless operation and easy charging.
- Steer and Brake Casters: Maneuver through tight hallways and lock the chair securely during procedures.
- Fold-Down Push Bar: Ergonomic handles provide caregivers with total control during

## SPECIFICATIONS

- Weight Capacity: 350 lbs
- Seat Height Range: 25" to 33"
- Backrest Range: 90° (upright) to 0° (flat)
- Seat Width: 24"
- Overall Width: 28.5" (rails down) to 31" (rails up)
- Backrest Width: 15" (Narrow for imaging)
- Standard Wheelbase: 30" x 20"
- Frame Material: Heavy-duty powder-coated steel

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)