

Champion Phlebotomy Chair | The Spirit



Stock Photo - Options May Vary

FACTORY NEW

The Champion Spirit Phlebotomy Chair is a high-performance medical chair that combines a modern aesthetic with rugged durability for daily clinical use. It provides a stable and secure platform for patients during lab work or vaccinations. You will notice how the sleek design fits perfectly into contemporary healthcare environments.

Consequently, your facility maintains a professional look while improving patient throughput. The convenient flip-up armrest, built into the Spirit Phlebotomy Chair, simplifies the entry and exit process for every patient. You can adjust the arm height to ensure the perfect angle for successful venipuncture.

This flexibility reduces physical strain for the phlebotomist and increases comfort for the patient.

FEATURES

- **Flip-Up Armrest:** Rotates out of the way to allow easy patient access from the side
- **Resin Construction:** Features a one-piece seat and backrest for fast, effective infection control
- **Adjustable Arm Height:** Tailors the draw position to the specific height of each patient
- **Optional Storage Drawer:** Keeps needles, alcohol pads, and bandages organized and accessible
- **Leveling Feet:** Ensures a rock-solid, wobble-free foundation on any floor surface

SPECIFICATIONS

- **Weight Capacity:** 300 lbs. (136 kg)
- **Overall Height:** 36" (91.4 cm)
- **Seat Dimensions:** 18" Width x 15" Depth x 20.5" Height
- **Arm Height Range:** 26.5" - 32.5" (67.3 cm - 82.6 cm)
- **Overall Width:** 29" (Standard) or 35.25" (With Drawer)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)