

Champion Phlebotomy Chair | The Solace



Stock Photo - Options May Vary

FACTORY NEW

The Champion Solace Phlebotomy Chair is a premium workstation that blends patient comfort with exceptional clinical utility. It provides a welcoming environment for blood draws and laboratory procedures. You can enhance the professional atmosphere of your clinic while ensuring patient safety.

Consequently, both staff and patients enjoy a more relaxed and efficient experience. The plush upholstered seat on the Solace, along with the high backrest, offers superior support during medical appointments. Patients will feel more at ease thanks to the ergonomic design and soft surfaces.

Furthermore, the flip-up armrest allows for effortless entry and exit for all body types.

FEATURES

- Upholstered Comfort: Features a padded seat and backrest for an enhanced patient experience
- Flip-Up Armrest: Swings upward to provide clear access for patients entering from the side
- Adjustable Arm Height: Customizes the draw site height to accommodate different clinician needs
- Optional Storage Drawer: Provides a convenient space for tubes, needles, and prep materials
- Adjustable Leveling Feet: Stabilizes the chair on various floor types for a secure foundation

SPECIFICATIONS

- Weight Capacity: 300 lbs. (136 kg)
- Overall Height: 44.5" (113 cm)
- Seat Dimensions: 18.75" Width x 18.5" Depth x 22.5" Height
- Arm Height Range: 27.75" - 33.75" (70.5 cm - 85.7 cm)
- Overall Width: 29" (Standard) or 35.25" (With Drawer)

**Brochure and images are a general overview of the product. Actual product options may vary.*