

Champion T4 Procedure Chair



Stock Photo - Options May Vary

FACTORY NEW

The Champion T4 Multi-Purpose Procedure Chair is a versatile medical platform that functions as a chair, a stretcher, and an exam table. It supports a wide range of procedures from minor surgeries to emergency wound care. Furthermore, you can adjust the height and backrest instantly with the push-button pendant.

Consequently, your staff saves time while maintaining precise control during critical tasks. Patient safety remains the highest priority of the T4 Procedure Chair with a 500-pound weight capacity. The triple column base ensures a rock-solid foundation for every patient.

Furthermore, the swing and tuck side rails allow for easy lateral transfers.

FEATURES

- Radiolucent Back: Allows for clear X-ray and imaging results while the patient remains seated.
- Swing & Tuck Side Rails: Provides unobstructed access for patient transfers and clinical procedures.
- Triple Column Base: Offers superior stability and a broad height range for ergonomic working positions.
- Central-Locking Brakes: Engages all wheels with one pedal to ensure a secure, stationary platform.
- Push-Button Pendant: Gives caregivers easy electronic control over all chair movements and positions.

SPECIFICATIONS

- Weight Capacity: 500 lbs. (226.8 kg)
- Seat Height Range: 23.5" to 43.5" (Standard Models)
- Seat Width: 24" (T4-B, T4-XB, T4-LB) or 28" (T4-WB)
- Overall Width (Rails Up): 31" to 35" depending on the model
- Patient Surface Length: 77" (Standard) or 71" (With footrest)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)