

Champion T5 Procedure Chair



Stock Photo - Options May Vary

FACTORY NEW

The Champion T5 Surgical Procedure Chair is a versatile medical platform that functions as a chair, stretcher, and operating table for specialty procedures including eye surgery, plastic surgery, and dental work. The dual-action headrest provides excellent access to the patient's head and neck with adjustable tilt and raise functions. Side rails swing down and tuck away for easy patient transfers, while the triple-post base ensures stability at all heights.

The chair features a central-locking brake system that secures all four wheels simultaneously, preventing unwanted movement during delicate procedures. A quick-release lever allows the backrest to flatten instantly during emergencies. The removable 24V lithium-ion battery eliminates floor clutter and can be swapped as needed.

Smooth maneuverability is enabled by specialized electrostatic steering casters that prevent static shock during transport between rooms. A push-button pendant provides electronic control over height, backrest angle, and Trendelenburg positioning. High-quality materials withstand frequent clinical cleaning and repeated use.

FEATURES

- Dual Articulating Headrest: Provides precise head positioning for ophthalmic and oral surgeries
- Triple Column Base: Delivers a broad height range and exceptional stability for surgical tasks
- Manual Quick-Release Back: Drops the backrest instantly for rapid emergency access
- Swing & Tuck Side Rails: Tucks completely out of the way for unobstructed patient access
- Central-Locking Brakes: Secures the entire unit firmly in place with a single pedal
- Surgical Support Leg: Adds extra stability to the base during weight-sensitive procedures
- Electrostatic Steering Caster: Ensures safe, static-free transport and easy maneuverability
- Push-Button Pendant: Offers electronic control over height, backrest, and Trendelenburg angles

SPECIFICATIONS

- Weight Capacity: 500 lbs. (227 kg)
- Seat Height Range: 23.75" to 43.5" (Standard Models)
- Seat Width: 24" (T5-B, T5-XB, T5-XLB) or 28" (T5-XWB)
- Backrest Articulation: 0° to 85°
- Trendelenburg Position: 20° (with 5° Reverse Trendelenburg)
- Headrest Adjustment: 30° - 60° articulation range
- Patient Surface Length: 75" (190.5 cm)
- Power: Removable 24V Lithium-Ion battery

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com