

# Winco Drop Arm Medical / Recovery / Dialysis Care Cliner w/Steel Casters | 6571XL



Stock Photo - Options May Vary

## FACTORY NEW

Winco Model 6571 XL Drop Arm Care Cliner. All of our best features, plus the added versatility of drop arms. Leave it to Winco to combine the best clinical chair features in the industry with the convenience of dual drop arms!

Our drop-arm Care Cliners make assisting your patients easier. In addition, our side panels lift, making clean-up quick and easy! The Winco 6571XL Features: Drop Arm Care Cliner has arms that can be lowered flush with the seat to aid in patient transfer.

It has non-porous blow-molded fold-down side trays on both sides. These bariatric recliners have a 450 lb weight capacity.

## FEATURES

- Ideal for dialysis, blood collection, oncology, respiratory, critical care, and cardiac
- 450 lb weight capacity
- Three-position recliner
- Infinite adjustable back
- Four heavy-duty 5" steel casters, positive locking swivel, and wheel brakes

## SPECIFICATIONS

- Weight: 111 lbs.
- Height: 49"
- Width with Tray Down: 37.75"
- Seat Height: 21"
- Seat Width: 25"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*