

Winco Drop Arm Medical / Recovery / Dialysis Recliner w/Steel Casters | 6551

FACTORY NEW

Winco 6551 Drop Arm CareCliner. All of our best features, plus the added versatility of drop arms. Leave it to Winco to combine the best clinical chair features in the industry with the convenience of dual drop arms!

Our drop-arm care cliners make assisting your patients easier. In addition, our side panels lift, making clean-up quick and easy. Winco Drop Arm Medical Recliner Features: The 6551 Drop Arm CareCliner has arms that can be lowered flush with the seat to aid in patient transfer.

It has non-porous blow-molded fold-down side trays on both sides. Each chair in the CareCliner Series features a standard Trendelenburg release with a handle accessible from either side.



Stock Photo - Options May Vary

FEATURES

- Ideal for dialysis, blood collection, oncology, respiratory, critical care, and cardiac
- Three-position recliner
- Infinite adjustable back
- Four heavy-duty Steel 5" casters/positive locking swivel and wheel brakes
- Exclusive 'one-touch" action drops arm securely in place

SPECIFICATIONS

- Weight: 105 lbs.
- Height: 49"
- Width: 33.75"
- Width with Tray Down: 33.75 "
- Seat Height: 21"

**Brochure and images are a general overview of the product. Actual product options may vary.*