

# Winco Extra Large Blood Drawing Chair with Vinyl Seat | 2575

## FACTORY NEW

Winco 2575 - Extra Large Blood Drawing Chair with Vinyl Seat. A strong, sturdy bariatric chair with all the features and quality you expect! A steel-constructed frame provides strength and durability for bariatric use.



Stock Photo - Options May Vary

The padded seat and back give exceptional patient comfort. Winco's commitment to quality will provide you with years of trouble-free service. Additional Features of Winco's 2575:

Winco's 2575 has all the same quality features as the 2572 chair, yet with an extra-large seat and back and added weight capacity for larger patients.

This is the same chair as our 2574 but without the drawer/storage cabinet. All available colors apply to this model.

## FEATURES

- Padded seat and back
- High-strength steel frame
- Non-slip rubber feet
- A neutral gray powder coat provides a durable finish
- Extra-large seat and back

## SPECIFICATIONS

- Seat Height: 22.5"
- Seat Width: 24.25"
- Overall Height: 36.5"
- Overall Width: 34"
- Armrest Adjustable Height: 28.5"-38.5"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)