

Winco Extra Large Padded Blood Drawing Chair with Cabinet | 2574

FACTORY NEW

Winco 2574 - Extra Large Padded Blood Drawing Chair with Cabinet. A strong, sturdy bariatric chair with all the features and quality you expect! The steel-constructed frame provides strength and durability for bariatric use. The padded seat and back give exceptional patient comfort.

Winco's commitment to quality will provide you with years of trouble-free service. Winco's 2574 offers the same quality features as the 2572, with an extra-large seat and back. It also provides an added weight capacity for larger patients.

All available colors apply to this model. Note: Looking for a Blood Chair with different features or weight capacity? Click the following link to view our catalog of Blood Draw Chairs.



Stock Photo - Options May Vary

FEATURES

- The steel-constructed frame provides strength and durability for bariatric use
- Padded seat and back give exceptional patient comfort
- Winco's commitment to quality will provide you with years of trouble-free service
- Extra large size
- Flip-up, contoured urethane arms

SPECIFICATIONS

- Height: 36.5"
- Width with Tray Down: 41.5"
- Length: 25"
- Weight: 72 lbs
- Seat Height: 22.5"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)